


HWS Screw Thread Units

"SILPET"

"SILPET PRO"

- more security for
your hose connection



H ▶ W ▶ S ▶ Screw Thread Units
SILPET / SILPET PRO

H ▶ W ▶ S ▶ Screw Thread Units

HWS-laboratory connection "MAXI"

enables the fitting of tubes
with a diameter from 8 or
10 mm for the main used
screw thread dimensions GL 14

complete with screw thread application
"SILPET"

	Art.-Nr.
GL 14 / 8 mm	1 264 908
GL 14 / 10 mm	1 264 910



HWS-screw thread application "SILPET"



by the combination of PTFE
(to the reaction side) and silicone
(pressure element) fast
and save screwing of
rod and tubes on GL-threads,
PTFE (TFM) highly flexible, free of pores,
silicone, permanent elastic
aging- and ultra violet rays-resistant
working temperature -50...+260 °C

- the screw thread that
keeps close

HWS-screw thread application "SILPET PRO"

a further development of the approved
"SILPET"-line. The silicone area of the
optimized seal holds tight all inserted
parts (like PTFE-hoses)
even under extreme conditions
like vacuum or pressure.
The ADS-security cap with a
inlaid silicone band avoids the
unwanted releasing of the cap
from the thread

- the maximum of security for you



screw thread application „SILPET“

for	Ø	Art.-Nr.	
	mm	(insert)	(complete)
GL14	1,6	1 264 141	1 064 1401
GL14	2,0	1 264 142	1 064 1402
GL14	3,0	1 264 143	1 064 1403
GL14	3,2	1 264 148	1 064 1403/2
GL14	4,0	1 264 144	1 064 1404
GL14	6,0	1 264 146	1 064 1406
GL18	1,6	1 264 181	1 064 1801
GL18	2,0	1 264 182	1 064 1802
GL18	3,0	1 264 183	1 064 1803
GL18	3,2	1 264 185	1 064 1803/2
GL18	4,0	1 064 184	1 064 1804
GL18	6,0	1 264 186	1 064 1806
GL18	8,0	1 264 188	1 064 1808
GL18	10	1 264 189	1 064 1810
GL25	3,0	1 264 253	1 064 2503
GL25	4,0	1 264 254	1 064 2540
GL25	6,0	1 264 256	1 064 2560
GL25	8,0	1 264 258	1 064 2580
GL25	10	1 264 260	1 064 2510
GL25	12	1 264 262	1 064 2512
GL25	14	1 264 264	1 064 2514
GL32	8,0	1 264 328	1 064 3208
GL32	10	1 264 330	1 064 3210
GL32	12	1 264 332	1 064 3212
GL32	14	1 264 334	1 064 3214
GL32	16	1 264 336	1 064 3216
GL32	18	1 264 338	1 064 3218
GL32	20	1 264 340	1 064 3220
GL45	6,0	1 264 456	1 064 4506
GL45	8,0	1 264 458	1 064 4508
GL45	10	1 264 460	1 064 4510
GL45	12	1 264 462	1 064 4512
GL45	14	1 264 464	1 064 4514
GL45	16	1 264 466	1 064 4516
GL45	18	1 264 468	1 064 4518
GL45	20	1 264 470	1 064 4520

application „SILPET-PRO“

GL14	1,6	1 265 141	1 065 1401
GL14	2,0	1 265 142	1 065 1402
GL14	3,0	1 265 143	1 065 1403
GL14	3,2	1 265 148	1 065 1403/2
GL14	4,0	1 265 144	1 065 1404
GL14	6,0	1 265 146	1 065 1406
GL18	1,6	1 265 181	1 065 1801
GL18	2,0	1 265 182	1 065 1802
GL18	3,0	1 265 183	1 065 1803
GL18	3,2	1 265 185	1 065 1803/2
GL18	4,0	1 265 184	1 065 1804
GL18	6,0	1 265 186	1 065 1806
GL18	8,0	1 265 188	1 065 1808
GL18	10	1 265 189	1 065 1810

please contact us in case of requirement

High Technology For Low - Cost Prices